



**CTL AMERICA**  
rethink what's possible®

## TAURUS™

### MODULAR THORACOLUMBAR PEDICLE SCREW SYSTEM

The Taurus™ Modular Thoracolumbar Pedicle System allows for true intra-operative flexibility. Headbodies can be connected in-situ or prior to insertion to allow direct screw to screw distraction for improved visualization. The patented set screw design eliminates head splay and cross threading to better assist set screw seating and rod reduction. The combination of intra-operative flexibility and system strength provides powerful benefits to both patients and surgeons.



Where strength connects with intra-operative flexibility

Modularity: Provides intra-operative flexibility

Helical Flange: Eliminates cross threading and head splay

Tension Head: Holds its position at any angle

## TAURUS PEDICLE SCREW SIZES

NON-CANNULATED DIAMETER	CANNULATED DIAMETER	SCREW LENGTH (millimeters)											
		30	35	40	45	50	55	60	70	80	90	100	110
5.5mm	5.5mm												
6.5mm	6.5mm												
7.5mm	7.5mm												
8.5mm	8.5mm												

## SYSTEM FEATURES

**Modularity:** Improves decortication of bone near the screw neck and optimizes graft placement for fusion

Headbodies can be connected in-situ or prior to insertion to allow direct screw to screw distraction for improved visualization

**Polyaxial:** Provides 50 degree of conical angulation

**Color-coded:** Assists in quicker screw identification

**Rod Options:** 5.5mm rods available in curved, straight, titanium and cobalt chrome



## SCREW FEATURES

**Dual Lead:** Enhances insertion speed and increases pullout strength

**Cannulation:** Enables guided screw placement (non-cannulated available)

**Self-Tapping:** Drives easily into bone while keeping the pedicle intact



## HEADBODY FEATURES

**Headbody Options:** Standard or reduction headbody

**Helical Flange:** Eliminates cross threading and head splay

**Tension Head:** Holds the headbody position at any angle